

9853 Pacific Heights Blvd. Suite D. San Diego, CA 92121, USA Tel: 858-263-4982

Email: info@abeomics.com

32-13520: ProSpec 8-Band Marker (discontinued)

Alternative Name: Protein Marker, Protein Ladder, Protein Standard, Mw Marker, Mw ladder, Mw Standard

Description

Source: E.Coli

Sterile Filtered blue solution.

ProSpec 8-Protein marker containing 15. 20, 25, 35, 50, 70,100 & 120 kDa. The recombinant proteins are produced in E.Coli is and purified by proprietary chromatographic techniques.

Product Info

Amount: 0.25 ml / 0.5 ml

Content:

ProSpec 8-protein marker solution (0.5ug/ul/protein) contains 12mM Tris-HCl pH-6.8, 5%

glycerol, 0.4% SDS 2.88mM 2-mercaptoethanol & 0.02% bromophenol blue.

Storage condition:

Storage condition:

Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Available traces that guides.

of time. Avoid multiple freeze-thaw cycles.